RAPIDVET™-H

CANE/CANINE (DEA 1.1)

GATTO/FELINE

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Blood group determination assay in dogs and cats

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The internationally accepted canine blood group system is the DEA (Dog Erythrocyte Antigen). Among the 8 antigens identified on the surface of the canine erythrocytes, DEA 1.1 is the most antigenic, involved in transfusions and in neonatal isoerythrolysis.

Dogs generally do not possess significant levels of isoantibodies to incompatible blood groups; they generally tolerate well an initial incompatible transfusion, although it stimulates the rapid production of antibodies. As a result, the half life of the transfused cells is very short and a second incompatible transfusion causes a deadly reaction.

Prior to mating it is very important to determine the blood type of the female and male dog, because the antibodies that are carried by the colostrum may cause neonatal isoerythrolysis in an incompatible newborn puppy.

RAPIDVET-H allows DEA 1.1 positive or negative dogs to be identified quickly and precisely, minimizing the risk of an incompatible transfusion and neonatal isoerythrolysis.

THE TEST The RapidVet-H assay is based on the hemagglutination reaction that occurs when an erythrocyte which contains a DEA 1.1 antigen on its surface membrane interacts with a monoclonal antibody proven specific to DEA 1.1, which is lyophilised on the Test Card.

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

INSTRUCTIONS

- Draw a minimum of 0.4 ml of blood from the patient into a syringe containing EDTA, as an anticoagulant.
- Dispense 1 drop (50 µl) of diluent solution into the well marked as DEA 1.1 positive Control and Patient Test, to reconstitute the lyophilised material on the card.
- Dispense 1 drop (50 µl) of Positive control and patient sample into the respective wells and mix the material for 10 seconds, using a new stirrer supplied with the kit.
- Rotate the card to mix reagents in each well and read results after 2 minutes.

CANE DEA 1.1 POS. / DOG DEA 1.1 POS

CANE DEA 1.1 NEG. / DOG DEA 1.1 NEG.

International Journal Articles Discussing our Products

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One blood system consisting of two antigens expressed either alone or in combination has been described in cats: Type A, Type B and Type AB. These types are important in transfusion reactions and in neonatal isoerythrolysis.

Cats possess isoantibodies to incompatible blood groups. About one-third of Type A cats have naturally occurring, low titered, anti-B antibodies. On the other hand all Type B cats have naturally occurring, highly titered anti-A antibodies. Cats with type B erythrocytes exhibit a catastrophic systemic reaction when transfused with type A blood because of their natural high-titered anti-A antibodies.

Those cats with type A erythrocytes, when first transfused with type B blood can exhibit either a mild reaction or develop moderate titers of anti-B antibodies that will result in a serious reaction upon a subsequent incompatible transfusion. The half life of the transfused cells will also be quite short.

Prior to mating it is very important to determine the blood type of the female and male cat, because the antibodies that are carried by the colostrum may cause in neonatal isoerythrolysis in type A newborn kittens, resulting from mating a type B queen with a type A male. The percentage of cats carrying the A and B antigen is shown in table 1.

RAPIDVET-H allows the blood group of the cat to be identified quickly and precisely, minimizing the risk of the incompatible transfusion and neonatal isoerythrolysis.

THE TEST The RapidVet-H (Feline) assay is based on the hemagglutination reaction that occurs when an erythrocyte, which contains either the type A or type B antigen or both (type AB) on its surface membrane, interacts with antibodies specific to the particular antigen, which is lyophilised on the Test Card.

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BREED	Blood frequency (%)		Allele Frequency (A+B=1)		Proportion of mating at risk for N.I. (%)
	TIPO A	TIPO B	A	B	
Abyssinian	86	14	.63	.37	12
Birman	84	16	.60	.40	13
British SH	60	40	.37	.63	24
Cornish Rex	66	34	.42	.58	23
Devon Rex	59	41	.36	.64	24
Domestic SH	99	1	.90	.10	1
Himalayan	93	7	.74	.26	6
Japanese Bobtail	84	16	.60	.40	13
Maine Coon	98	2	.86	.14	2
Norwegian Forest	93	7	.74	.26	6
Persian	86	14	.63	.37	12
Scottish Fold	82	18	.58	.42	15
Siamese	100	0	1.00	.00	0
Sphinx	82	18	.56	.44	16
Somali	83	17	.59	.41	14
Tonkinese	100	0	1.00	.00	0

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

INSTRUCTIONS

- Draw a minimum of 0.4 ml of blood from the patient into a syringe containing EDTA, as an anticoagulant.
- Dispense 1 drop (50 µl) of diluent solution into the wells marked as Type A and Type B, to reconstitute the lyophilised material on the card.
- Dispense 1 drop (50 µl) of patient sample into the two respective wells and mix the material for 10 seconds, using a new stirrer supplied with the kit.
- Rotate the card to mix reagents in each well and read results after 2 minutes.

GATTO TIPO A / CAT TYPE A

GATTO TIPO B / CAT TYPE B

Free phone number for Italy:

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Internet: <http://www.agrolabo.it> - E-mail: marketing@agrolabo.it

Head Office, Laboratory and Production Centre: Agrolabo S.p.A., via Masero snc - 10010 Scarmagno (TO) - Italy - Tel. 0039 0125 731111 - Fax 0039 0125 731190